




Reading your results - Lateral Flow Test (LFT)

NAME:

DATE:

What did my test show?	What does my result mean?
<p style="text-align: center;">Negative</p>  <p style="text-align: center;">Only the control 'C' line appeared in the results window</p>	<p>A visible control 'C' line means your test device worked properly</p> <p>The absence of the 'T' line indicates that the test is negative</p> <p>A negative result means that it is less likely that you have COVID -19</p> <p>Safety measures remain important because this test is not as sensitive as a PCR test. You should continue to adhere to all current public health guidance</p> <p>If you develop symptoms you should self-isolate immediately and request a PCR test in line with current public health guidance by calling the Coronavirus helpline on 0800 735 5566.</p>
<p style="text-align: center;">Positive</p>  <p style="text-align: center;">The T and control 'C' lines appeared in the results window</p>	<p>A visible control 'C' line means your test device worked properly</p> <p>A visible 'T' line indicates that the test is positive</p> <p>You will need to go home and begin self-isolation immediately. You should isolate within the home away from other household members. If you are unable to isolate within the home away from other household members, your whole household will need to isolate until you receive a negative PCR result.</p> <p>Call the Coronavirus helpline to arrange a home visit for a PCR test.</p> <p>If the PCR test is negative, you can return to work and your household will no longer need to isolate</p> <p>If your PCR test is negative you can re-join the Covid-19 Community Testing Programme (Lateral Flow Testing) immediately.</p> <p>If your PCR test is positive you must wait 28 days from the last day of isolation before re-joining the Covid-19 Community Testing programme.</p>
<p style="text-align: center;">Invalid or Void</p>  <p style="text-align: center;">The control 'C' line is absent from the results window</p>	<p>The absence of a 'C' line means that your test device did not work properly - your test is invalid or void</p> <p>You can continue working as normal but must arrange another LFT at the earliest possible opportunity</p>

Please note: LFT provides an additional layer of precaution against Covid-19 outbreaks which supports existing public health measures. A negative LFT result does NOT categorically indicate that you are COVID-19 negative. Please continue to strictly adhere to public health guidance even when presented with a negative result.

If you have any questions regarding your result, please contact the Coronavirus helpline on **0800 735 5566**.